

Items classified as SHARPS WASTE:
Needles • Syringes • Blades • Scalpels
• Clinical Glass or any objects that can cause a cut, abrasion or puncture wound.

Do's:

- ✓ Ensure the lid is secured properly along the perimeter of the container
- ✓ Place sharps containers in the mounted bracket or secured area
- ✓ Ensure that only sharp items are placed in the sharps container (refer to above classification)
- ✓ Securely lock the flaps on top of the sealed lid once the sharps reaches the 'fill line'
- ✓ Once filled, place sealed sharps in designated waste holding / storage area
- ✓ Always keep sharps containers in an upright position

Don'ts:

- ✗ Never exceed the fill line indicated on the sharps container
- ✗ Do not leave sharps containers in areas where unauthorized persons may have access to it
- ✗ Do not dispose of paper, gloves and general non-sharp HCRW into sharps containers

Please note that BCL cannot remove any sharps containers that have not been sealed and secured correctly.

* Note:

All containers are not proportioned to each other



20 Litre



7.6 Litre



6.6 Litre



5 Litre



1 Litre

SEPARATION OF CONTAMINATED MEDICAL WASTE AT SOURCE IS ESSENTIAL TO ENSURE THAT THE SPECIFIC WASTE STREAMS ARE PLACED INTO THE CORRECT CONTAINERS. DO NOT DISPOSE OF ANY GENERAL WASTE INTO MEDICAL WASTE CONTAINERS.